



Quick Start Guide

RxKeeper®
Designed in USA
Platform Maintained in USA
USA FDA Registered
USA Patent Issued

Let's Get Started

RxKeeper® is a sturdy and secure container designed to hold your medications and vitamins. It has a temperature display and a motion-detecting alarm clock that will ring and flash if you forget to take your medication.

Open the RxKeeper® Container, take out the keys, and store the keys in a safe place.

Please DO NOT discard the Red RxKeeper® cardboard box.

Store your medications and supplements inside your RxKeeper® Container and take your medication at the prescribed time.



Note: After locking the container, please keep one key in a safe place.

Setting the RxKeeper® Alarm Clock

The RxKeeper® Alarm Clock (shown below) has three red buttons.



To turn ON, the AM alarm PRESS AM Alarm for 3 seconds.

To turn ON, the PM alarm PRESS PM Alarm for 3 seconds.

Instructions for Setting the RxKeeper® Alarm Clock

Action 1 Setting Device Time	Action 2 Setting AM Alarm	Action 3 Setting PM Alarm
1. After powering the Alarm Clock, press and hold Time for 3 seconds. The "hour" digits of the device time will flash.	1. Press Time . The "AM Alarm Hour" digits will flash.	1. Press Time . The "PM Alarm Hour" digits will flash.
2. Press PM once to increase the hour digits, or press AM once to decrease the hour digits. Note: Hold for 3 seconds to increase or decrease the hour digits at a faster rate.	2. Press PM once to increase the hour digits, or press AM once to decrease the hour digits. Note: Hold for 3 seconds to increase or decrease the hour digits at a faster rate.	2. Press PM once to increase the hour digits, or press AM once to decrease the hour digits. Note: Hold for 3 seconds to increase or decrease the hour digits at a faster rate.
3. Set the hour digits.	3. Set the AM Alarm hour digits.	3. Set the PM Alarm hour digits.
4. Press Time . The "minute" digits of the device time will flash.	4. Press Time . The "AM Alarm Minute" digits will flash.	4. Press Time . The "PM Alarm Minute" digits will flash.
5. Press PM once to increase the minute digits, or press AM once to decrease the minute digits. Note: Hold for 3 seconds to increase or decrease the minute digits at a faster rate.	5. Press PM once to increase the minute digits, or press AM once to decrease the minute digits. Note: Hold for 3 seconds to increase or decrease the minute digits at a faster rate.	5. Press PM once to increase the minute digits, or press AM once to decrease the minute digits. Note: Hold for 3 seconds to increase or decrease the minute digits at a faster rate.
6. Set the minute digits for the device.	6. Set the AM Alarm Minute digits.	6. Set the PM Alarm Minute digits.
Note: You must perform Action 1,2,3 sequentially. Pressing the Time button after each Time setting to sequence through device time, AM alarm time, and PM alarm time.		7. After all these steps, click TIME to confirm all settings.

Replacing the RxKeeper® Device Battery

ONLY when a Care coordinator asks you to change the battery:



1. Locate the "Sensor Device" inside the RxKeeper® container.
2. Pull the device up and detach it slowly from the container tracks using both hands. (Fig. 1)
3. Remove the black plug. (Fig. 2)
4. Pull the clamps out and remove the black back cover.

Insert two new AA batteries. (Fig 3)



5. Put the clamps back. Insert the plug back securely.
6. Insert the "Sensor Device" back into the RxKeeper container tracks.



Safety Information

The RxKeeper® platform and medication container are intended as a backup medication management system that can be customized for each user. It is not the primary tool for determining whether and when to take medication. You must consult with healthcare providers and make medical decisions based on their advice. This platform may not guarantee the on-time administration of medication or supplement alerts. The user is responsible for consulting professional medical advisors. We do not accept any liability concerning any activities you may undertake using the RxKeeper® platform services. Status Alert LLC will not be liable for any injury resulting or arising from, or related to, the use of the services (even if unauthorized) whatsoever.

IF YOU RELY ON ANY CONTENT FOUND ON THE DEVICE OR SERVICES, YOU DO SO SOLELY AT YOUR OWN RISK. YOUR USE OF THE SERVICES AND DEVICE IS NOT A SUBSTITUTE FOR HEALTHCARE. NO DOCTOR-PATIENT OR PROVIDER-PATIENT RELATIONSHIP IS CREATED WHEN YOU USE THE SERVICES OR CONTENT. THIS IS TRUE WHETHER SUCH CONTENT IS PROVIDED BY OR THROUGH THE USE OF THE SERVICES OR THROUGH ANY OTHER COMMUNICATIONS FROM STATUS ALERT LLC.



WARNING: Keep RxKeeper® out of the reach of children.

WARNING: Keep RxKeeper® Container keys in a safe place.

CAUTION: Call your doctor for any medical issues.

CAUTION: Call 911 in case of an Emergency.